



Product Name: Recombinant Human CDHR5 Protein
Catalog #: PHH4086



Summary

Name	CDHR5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q9HBB8
Host	HEK293
Species	Human
Molecular Mass	110 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SVNKDIFEVE ENTNVTEPLV DIHVPEGQEV TLGALSTPFA FRIQGNQLFL NVTPDYEEKS LLEAQLLCQS GGTLVTQLRV FVSVLDVNDN APEFPFKTKE IRVEEDTKVN STVIPETQLQ AEDRDKDDIL FYTLQEMTAG ASDYFSLVSV NRPALRLDRP LDFYERPNMT FWLLVRDTPG ENVEPSHTAT ATLVLVNVVPA DLRPPWFLPC TFSDBGYVCIQ AQYHGAVPTG HILPSPLVLR PGPIYAEDGD RGINQPIYS IFRGNVNGTF IHPDSGNLT VARSVPSPMT FLLLVKGQQA DLARYSVTQV TVEAVAAAGS PPRFPQRLYR GTVARGAGAG VVVKDAAAPS QPLRIQAQDP EFDLNSAIT YRITNHSHFR MEGEVVLTST TLAQAGAFYA EVEAHNTVTS GTATTVIEIQ VSEQEPPSTD VPPSPEAGGT TGPWTSTTSE VPRPPEPSQG PSTTSSGGGT GPHPPSGTTL RPPTSSTPGG PPGAENSTSH QPATPGGDTA QTPKPGTSQP MPPGVGTSTS HQPATPSGGT AQTPPEPGTSQ PMPPSMGTST SHQPATPGGG TAQTPEAGTS QPMPPGMGTS TSHQPTTPGG GTAQTPEPGT SQPMPLSKST PSSGGGPPSED

image



Product Name: Recombinant Human CDHR5 Protein
Catalog #: PHH4086



Background

Alternative Names

Cadherin-related family member 5 (Mu-protocadherin) (Mucin and cadherin-like protein) (Mucin-like protocadherin) (MLPCDH)

Note

For Research Use Only , Not for Diagnostic Use