



Product Name: Recombinant Human CDHR5 Protein
Catalog #: PHH4087



Summary

Name	CDHR5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q9HBB8
Host	HEK293
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VPRPPEPSQG PSTTSSGGGT GPHPPSGTTL RPPTSSTPGG PPGAENSTSH QPATPGGDTA QTPKPGTSQP MPPGVGTSTS HQPATPSGGT AQTPEPGTSQ

image



Background

Alternative Names Cadherin-related family member 5 (Mu-protocadherin) (Mucin and cadherin-like protein) (Mucin-like protocadherin) (MLPCDH)

Note

For Research Use Only , Not for Diagnostic Use