



**Product Name: Recombinant Human NECTIN2 Protein**  
**Catalog #: PHH4110**



## Summary

<b>Name</b>	NECT2_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q92692
<b>Host</b>	HEK293
<b>Species</b>	Human
<b>Molecular Mass</b>	36 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	<p>EVRGQLGGTV ELPCHLLPPV PGLYISLVTW QRPDAPANHQ NVAAFHPKMG PSFSPKPGS ERLSFVSAKQ STGQDTEAEL QDATLALHGL TVEDEGNYTC EFATFPKGSV RGMTWLRVIA KPKNQAEAQK VTFSQDPTTV ALCISKEGRP PARISWLSSL DWEAKETQVS GTLAGTVTVT SRFTLVPSGR ADGVTVTCKV EHESFEEPAL IPVTLVRYYP PEVSISGYDD NWYLGRTDAT LSCDVRSNPE PTGYDWSTTS GTFPTSAVAQ GSQLVIHAVD SLFNNTTFVCT VTNAVGMGRA EQVIFVRETP NTAGAGATGG</p>

## image



## Background

**Alternative Names** Nectin-2 (Herpes virus entry mediator B) (Herpesvirus entry mediator B) (HveB) (Nectin cell adhesion molecule 2) (Poliovirus receptor-related protein)

**Product Name: Recombinant Human NECTIN2 Protein**  
**Catalog #: PHH4110**



---

2) (CD antigen CD112)

**Note**

For Research Use Only , Not for Diagnostic Use