


Product Name: Recombinant Human VSIG4 Protein
Catalog #: PHH4116



Summary

Name	VSIG4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y279
Host	HEK293
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PILEVPESVT GPWKGDVNLP CTYDPLQGYT QVLVKWLVR GSDPVTIFLR DSSGDHIQQA KYQGRLHVSH KVPGDVSLQL STLEMDDRSH YTCEVTWQTP DGNQVVRDKI TELRVQKLSV SKPTVTTGSG YGFTVPQGMR ISLQCQARGS PPISYIWYKQ QTNNQEPIKV ATLSTLLFKP AVIADSGSYF CTAKGQVGSE QHSDIVKFVV KDSSKLLKTK TEAPTTMTYP LKATSTVKQS WDWTTDMDGY LGETSAGPGK

image



Background

Alternative Names	V-set and immunoglobulin domain-containing protein 4 (Protein Z39Ig)
--------------------------	--

Product Name: Recombinant Human VSIG4 Protein
Catalog #: PHH4116



Note

For Research Use Only , Not for Diagnostic Use