


Product Name: Recombinant Human XG Protein
Catalog #: PHH4126



Summary

Name	XG_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	P55808
Host	HEK293
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LDDPEPTKKP NSDIYKPKPK PYYPQENPD SGGNIYPRPK PRPQPQPGNS GNSGGYFNDV DRDDGRYPPR

image



Background

Alternative Names Glycoprotein Xg (Protein PBDX)

Note

For Research Use Only , Not for Diagnostic Use