


Product Name: Recombinant Human ICAM5 Protein
Catalog #: PHH4130



Summary

Name	ICAM5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UMF0
Host	HEK293
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VGSERPVSCT LDGLFPASEA RVYLALGDQN LSPDVTLEGD AFVATATATA SAEQEGARQL VCNVTLGGEN RETRENTIY SFPAPLLTLS EPSVSEGMV TVTCAAGAQA LVTLEGVPAA VPGQPAQLQL NATENDRRS FFCDATLDVD GETLIKNRSA

image



Background

Alternative Names Intercellular adhesion molecule 5 (ICAM-5) (Telencephalin)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ICAM5 Protein
Catalog #: PHH4130

