


**Product Name: Recombinant Human CD84 Protein**  
**Catalog #: PHH4132**



## Summary

<b>Name</b>	SLAF5_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q9UIB8
<b>Host</b>	HEK293
<b>Species</b>	Human
<b>Molecular Mass</b>	24 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	GILGESVTFP VNIQEPRQVK IIAWTSKTSV AYVTPGDSET APVVTVTHRN YYERIHALGP NYNLVISDLR MEDAGDYKAD INTQADPYTT TKRYNLQIYR RLGKPKITQS LMASVNSTCN VTLTCSVEKE EKNVTYNWSP LGEEGNVLQI FQTPEDQELT YTCTAQNPVS NNSDSISARQ LCADIAMGFR

## image



## Background

<b>Alternative Names</b>	SLAM family member 5 (Cell surface antigen MAX.3) (Hly9-beta) (Leukocyte differentiation antigen CD84) (Signaling lymphocytic activation molecule 5) (CD antigen CD84)
--------------------------	--

## Note

**Product Name: Recombinant Human CD84 Protein**  
**Catalog #: PHH4132**



---

For Research Use Only , Not for Diagnostic Use