



Product Name: Recombinant Human CD81 Protein
Catalog #: PHH4135



Summary

Name	CD81_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-His & N-hFc
Accession #	P60033
Host	HEK293
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KDVKQFYDQA LQQAVVDDDA NNAKAVVKTF HETLDCCGSS TLTALTTSVL KNNLCPSGSN IISNLFKEDC

image



Background

Alternative Names	CD81 antigen (26 kDa cell surface protein TAPA-1) (Target of the antiproliferative antibody 1) (Tetraspanin-28) (Tspan-28) (CD antigen CD81)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use