


Product Name: Recombinant Human NRP2 Protein
Catalog #: PHH4143



Summary

Name	NRP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	O60462
Host	HEK293
Species	Human
Molecular Mass	55 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VETLGPTVKS EETTPYPTE EEATECGENC SFEDDKDLQL PSGFNCNFD LEPCGWMYD HAKWLRTTWA SSSPNDRTF PDDRNFLRLQ SDSQREGQYA RLISPPVHLP RSPVCMFQY QATGGRGVAL QVVREASQES KLLWVIREQ GGEWKHGRII LPSYDMEYQI VFEGVIGKGR SGEIAIDDIR ISTDVPLENC MEPISAFAGE NFKVDIPEIH EREGYEDEID

image



Background

Alternative Names	Neuropilin-2 (Vascular endothelial cell growth factor 165 receptor 2)
--------------------------	---

Note

Product Name: Recombinant Human NRP2 Protein
Catalog #: PHH4143



For Research Use Only , Not for Diagnostic Use