



## Summary

<b>Name</b>	IL4RA_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P24394
<b>Host</b>	HEK293
<b>Species</b>	Human
<b>Molecular Mass</b>	25 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	EPTCVSDYMS ISTCEWKMG PTNCSTELRL LYQLVLLSE AHTCIPENNG GAGCVCHLLM DDVVSADNYT LDLWAGQQLL WKGSEFKPSEH VKPRAPGNLT VHTNVSDTLL LTWSNPYPPD NYLYNHLYA VNIWSENDPA DFRIYNVTYL EPLRIAASST LKSGISYRAR VRAWAQCYN T WSEWSPSTK WHNSYREPFE

## image



## Background

<b>Alternative Names</b>	Interleukin-4 receptor subunit alpha (IL-4 receptor subunit alpha) (IL-4R subunit alpha) (IL-4R-alpha) (IL-4RA) (CD antigen CD124) [Cleaved into: Soluble interleukin-4 receptor subunit alpha (Soluble IL-4 receptor subunit alpha) (Soluble IL-4R-alpha) (sIL4Ralpha/prot) (IL-4-binding protein) (IL4-BP)]
--------------------------	---

**Product Name: Recombinant Human IL4R Protein**  
**Catalog #: PHH4146**



---

**Note**

For Research Use Only , Not for Diagnostic Use