



Product Name: Recombinant Human SHH Protein
Catalog #: PHH4168



Summary

Name	SHH_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-His & N-MBP & N-His
Accession #	Q15465
Host	HEK293
Species	Human
Molecular Mass	63 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GKRRHPKKLT PLAYKQFIPN VAEKTLGASG RYEGKISRNS ERFKELTPNY NPDIIFKDEE NTGADRLMTQ RCKDKLNALA ISVMNQWPGV KLRVTEGWDE DGHHSSESLH YEGRAVDITT SDRDRSKYGM LARLAVEAGF DWVYYESKAH IHCSVKAENS

image



Background

Alternative Names	Sonic hedgehog protein (SHH) (HHG-1) (Shh unprocessed N-terminal signaling and C-terminal autoprocessing domains) (ShhNC) [Cleaved into: Sonic hedgehog protein N-product (ShhN) (Shh N-terminal processed signaling domains) (ShhNp)]
--------------------------	--

Product Name: Recombinant Human SHH Protein
Catalog #: PHH4168



Note

For Research Use Only , Not for Diagnostic Use