


Product Name: Recombinant Human MMP15 Protein
Catalog #: PHH4177



Summary

Name	MMP15_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P51511
Host	HEK293
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QPTQLPTVT PRRPGRPDHR PRRPPQPPPP GGKPERPPKP GPPVQPRATE RPDQYGPNIC DGDFDTVAML

image



Background

Alternative Names	Matrix metalloproteinase-15 (MMP-15) (EC 3.4.24.-) (Membrane-type matrix metalloproteinase 2) (MT-MMP 2) (MTMMP2) (Membrane-type-2 matrix metalloproteinase) (MT2-MMP) (MT2MMP) (SMCP-2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use