


Product Name: Recombinant Human TIMD4 Protein
Catalog #: PHH4196



Summary

Name	TIMD4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q96H15
Host	HEK293
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VLPSDSWSSV ESTSADTVLL TSKESKVWDL PSTSHVSMWK TSDSVSSPQP GASDTAVPEQ NKTTKTGQMD

image



Background

Alternative Names T-cell immunoglobulin and mucin domain-containing protein 4 (TIMD-4) (T-cell immunoglobulin mucin receptor 4) (TIM-4) (T-cell membrane protein 4)

Note

For Research Use Only , Not for Diagnostic Use