


## Summary

<b>Name</b>	CO6A2_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P12110
<b>Host</b>	HEK293
<b>Species</b>	Human
<b>Molecular Mass</b>	107 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	QQQEISPDT TERNNNCPEK TDCPIHVYFV LDTSESVTMQ SPTDILLFHM KQFVPQFISQ LQNEFYLDQV ALSWRYGGLH FSDQVEVFSP PGSDRASFIK NLQGISSFRR GTFTDCALAN MTEQIRQDRS KGTVHFAVVI TDGHVTGSPC GGIKLQAERA REEGIRLFAV APNQLKEQG LRDIASPHE LYRNDYATML PDSTEIDQDT INRIIKVMKH EAYGECYKVS CLEIPGPSGP KGYRGQKGAK GNMGEPGEPG QKGRQGDPGI EGPIGFPGPK GVPGFKGEKG EFGADGRKGA PGLAGKNGTD GQKGLGRIG PPGCKGDPGN RGPDPYPGEA GSPGERGDQG GKGDPGRPGR RGPPGEIGAK GSKGYQGNSG APGSPGVKGA KGGPGPRGPK GEPGRRGDGP TKGSPGSDGP KGEKGDGPPE GPRGLAGEVG NKGAKGDRGL PGPRGPQ GAL GEPGKQGSRG DPGDAGPRGD SGQP GPKGDP GRPGFSYPGP RGAPGEKGEP GPRPEGGRG DFGLKGEPGR KGEKGEPADP GPPGEPGPRG PRGVPGPEGE PGPPGDPGLT ECDVMTYVRE TCGCCDCEKR CGALDVVFI DSSESIGYTN FTLEKNFVIN VVNRLGAIK DPKSETGTRV GVVQYSHEGT FEAIQLDDER IDSLSSFKEA VKNLEWIAGG TWTPSALKFA YDRLIKESRR QKTRVFAVVI TDGRHDPRDD DLNLRALCDR DVTVTAIGIG DMFHEKHESE NLYSIACDKP QQVRNMTLFS DLVAEKFIDD MEDVLCPPDQ IVCPLPCQT ELSVAQCTQR PVDIVFLLDG SERLGEQNFH KARRFVEQVA RRLTLARRDD DPLNARVALL QFGGPGEQQV AFPLSHNLTA IHEALETQY LNSFSHVAG VVHAINAIVR SPRGGARRHA ELSFVFLTDG VTGNDLHES AHSMRKNVV PTVLALGSDV DMDVLTTLISL GDRAAVFHEK DYDSLAPQPGF FDRFIRWIC

## image

**Product Name: Recombinant Human COL6A2 Protein**  
**Catalog #: PHH4201**



## Background

**Alternative Names**

Collagen alpha-2(VI) chain

## Note

For Research Use Only , Not for Diagnostic Use