



Product Name: Recombinant Human CD164 Protein
Catalog #: PHH4214



Summary

Name	MUC24_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q04900
Host	HEK293
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FWIECKDESY CSHNSTVSDC QVGNTTDFCS VSTATPVPTA NSTAKPTVQP SPSTTSKTVT TSGTTNNTVT PTSQPVRKST

image



Background

Alternative Names	Sialomucin core protein 24 (MUC-24) (Endolyn) (Multi-glycosylated core protein 24) (MGC-24) (MGC-24v) (CD antigen CD164)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use