


Summary

Name	ECSCR_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-GST & C-His
Accession #	Q19T08
Host	HEK293
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TSSSQGGLGG LSLTTEPVSS NPGYIPSEA NRPSHLSSTG TPGAGVPSSG RDGGTSRDTF QTVPNSTTM SLSMREDATI LPSPTSETVL

image



Background

Alternative Names Endothelial cell-specific chemotaxis regulator (Apoptosis regulator through modulating IAP expression) (ARIA) (Endothelial cell-specific molecule 2)

Note

For Research Use Only , Not for Diagnostic Use