



Product Name: Recombinant Human MEGF9 Protein
Catalog #: PHH4223



Summary

Name	MEGF9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q9H1U4
Host	HEK293
Species	Human
Molecular Mass	42 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GEP SHPFRA TAPTAQAPRT GPPRATVHRP LAATSPAQSP ETTPLWATAG PSSTTFQAPL GPSPTTPPAA ERTSTTSQAP TRPAPTTLST TTGPAPTTTPV ATTVPAPTTTP RTPTPDLPS SNSSVLPTPP

image



Background

Alternative Names	Multiple epidermal growth factor-like domains protein 9 (Multiple EGF-like domains protein 9) (Epidermal growth factor-like protein 5) (EGF-like protein 5)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MEGF9 Protein
Catalog #: PHH4223

