

Product Name: Recombinant Human CPXM1 Protein
Catalog #: PHH4240



Summary

Name	CPXM1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96SM3
Host	HEK293
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QPGTTKVPGS TPALHSSPAQ PPAETANGTS EQHVRIRVIK KKKVIMKKRK KLTLTRPTPL VTAGPLVTPT PAGTLDPAEK QETGCPPLGL ESLRVSDSRL EASSQSFGGL GPHRGRLNIQ SGLEDGDLYD GAWCAEEQDA DPWFQVDAGH PTRFSGVITQ GRNSVWRYDW VTSYKVQFSN DSRTWWGSRN HSSGMDAVFP ANSDPETPVL NLLPEPQVAR FIRLLPQTWL QGGAPCLRAE

image



Background

Alternative Names Probable carboxypeptidase X1 (EC 3.4.17.-) (Metalloproteinase CPX-1)

Note

Product Name: Recombinant Human CPXM1 Protein
Catalog #: PHH4240



For Research Use Only , Not for Diagnostic Use