


## Summary

<b>Name</b>	SEM6C_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q9H3T2
<b>Host</b>	HEK293
<b>Species</b>	Human
<b>Molecular Mass</b>	63 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	LPLLISDLQG TSPLSWFRGL EDDAVAAELG LDFQRFLTIN RTLLVAARDH VFSFDLQAE EGEGLVPNKY LTWRSQDVEN CAVRGKLTDE CYNIRVLP WDSQTLACG TNSFSPVCRS YGITSLQQEG EELSGQARCP FDAQSNVAI FAEGSLYSAT AADFQASDAV VYRSLGPQPP LRSKYDSKW LREPHVQAL EHGDHVYFFF REVSVEDARL GRVQFSRVAR VCKRDMGGSP RALDRHWTSF LKLRLNCSVP GDSTFYFDVL QALTGPVNLH GRSALFGVFT TQNSIPGSA VCAFYLDEIE RGFEGKFKEQ RSLDGAWTPV SEDRVPSRP GSCAGVGGAA LFSSRDLPD DVLTFIKAHP LLDPAVPPVT HQPLLTLSR ALLTQVAVDG MAGPHSNITV MFLGSNDGTV LKVLTPGGRS GGPEPILLEE IDAYSPARCS GKRTAQTARR IIGLELDTEG HRLFVAFSGC IVYPLSRCA RHGACQRSCL ASQDPYCGWH SSRGCVDIRG SGGTDVDQAG NQESMEHGDC QDGATGSQSG PGDSAYGVRR DLPPASASRS

## image



**Product Name: Recombinant Human SEMA6C Protein**  
**Catalog #: PHH4246**



---

## Background

**Alternative Names**                      Semaphorin-6C (Semaphorin-Y) (Sema Y)

## Note

For Research Use Only , Not for Diagnostic Use