


Product Name: Recombinant Human MICA Protein
Catalog #: PHH4257



Summary

Name	MICA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q29983
Host	HEK293
Species	Human
Molecular Mass	75 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NLTVLSWDGS VQSGFLTEVH LDGQPFLRCD RQKCRAPQG QWAEDVLGNK TWDRETRDLT GNGKDLRMTL AHIKDQKEGL HSLQEIRVCE IHEDNSTRSS QHFYYDGELF LSQNLETKEW TMPQSSRAQT LAMNVRNFLK EDAMKTKTHY HAMHADCLQE LRRYLKSGVV LRRTVPPMVN VTRSEASEGN ITVTCRASGF YPWNITLSWR QDGVSLSHDT QQWGDVLPDG NGTYQTWVAT RICQGEEQRF TCYMEHSGNH STHPVPSGKV

image



Background

Alternative Names	MHC class I polypeptide-related sequence A (MIC-A)
--------------------------	--

Product Name: Recombinant Human MICA Protein
Catalog #: PHH4257



Note

For Research Use Only , Not for Diagnostic Use