


Product Name: Recombinant Human SLITRK1 Protein
Catalog #: PHH4279



Summary

Name	SLIK1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q96PX8
Host	HEK293
Species	Human
Molecular Mass	52 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FLRDNKHISI RKSHFVDYKN LILLDLGNNN IATVENNTFK NLLDLRWLYM DSNYLDTLSR EKFAQLQNLE

image



Background

Alternative Names SLIT and NTRK-like protein 1 (Leucine-rich repeat-containing protein 12)

Note

For Research Use Only , Not for Diagnostic Use