



Product Name: Recombinant Human PZP Protein
Catalog #: PHH4302



Summary

Name	PZP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-GST & C-His
Accession #	P20742
Host	HEK293
Species	Human
Molecular Mass	44 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SEKFEIENCL ANKVDLSFSP AQSPASHAH LQVAAAPQSL CALRAVDQSV LLMKPEAELS VSSVYNLLTV KDLTNFPDNV DQQEEEQGHC PRPFFIHNGA IYVPLSSNEA DIYSFLKGMG LKVFTNSKIR KPKSCSVIPS VSAGAVGQGY

image



Background

Alternative Names Pregnancy zone protein (C3 and PZP-like alpha-2-macroglobulin domain-containing protein 6)

Note

For Research Use Only , Not for Diagnostic Use