



Product Name: Recombinant Human RARRES1 Protein
Catalog #: PHH4304



Summary

Name	TIG1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P49788
Host	HEK293
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QAARAALHFF NFRSGSPSAL RVLAEVQEGR AWINPKEGCK VHVVFSTERY NPESLLQEGE GRLGKCSARV FFKNQKPRPT INVTCTRLIE KKKRQQEDYL LYKQMKQLKN PLEIVSIPDN HGHIDPSLRL IWDLAFLGSS YVMWEMTTQV SHYYLAQLTS VRQWKTNDDT IDFDYTVLLH ELSTQEIIPC RIHLVWYPGK PLKVKYHCQE LQTPEEASGT EEGSAVVPTE

image



Background

Alternative Names	Retinoic acid receptor responder protein 1 (Phorbol ester-induced gene 1 protein) (PERG-1) (RAR-responsive protein TIG1) (Tazarotene-induced gene 1 protein)
--------------------------	--

Product Name: Recombinant Human RARRES1 Protein
Catalog #: PHH4304



Note

For Research Use Only , Not for Diagnostic Use