



## Summary

<b>Name</b>	APCD1_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	N-His
<b>Accession #</b>	Q8J025
<b>Host</b>	HEK293
<b>Species</b>	Human
<b>Molecular Mass</b>	55 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	DSRSHPRSLE KSAWRAFKES QCHHMLKHLH NGARITVQMP PTIEGHWVST GCEVRSQPEF ITRSYRFYHN NTFKAYQFY GSNRCTNPTY TLIIRGKIRL RQASWIIRGG TEADYQLHNV QVICHTEAVA EKLQQVNRT CPGFLADGGP WVQDVAYDLW REENGCECTK AVNFAMHELQ LIRVEKQYLH HNL DHLVEEL FLGDIHTDAT QRMFYRPSSY QPPLQNAKNH DHACIACRII YRSDEHHPPI LPPKADLTIG LHGEWVSQRC EVRPEVFLT RHFIFHDNNN TWEGHYHYHS DPVCKHPTFS IYARGRYSRG VLSSRVMGGT EFVFKVNHMK VTPMDAATAS LLNVFNGNEC GAEGSWQVGI QQDVHTHTNGC VALGIKLPHY EYEIFKMEQD ARGRYLLFNG QRPSDGSSPD RPEKRATSYQ MPLVQCASSS PRAEDLAEDS GSSLYGRAPG

## image



## Background

**Product Name: Recombinant Human APCDD1 Protein**  
**Catalog #: PHH4324**



---

**Alternative Names**

Protein APCDD1 (Adenomatosis polyposis coli down-regulated 1 protein)

**Note**

For Research Use Only , Not for Diagnostic Use