


Summary

Name	DHH_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43323
Host	HEK293
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GRRRYARKQL VPLLYKQFVP GVPERTLGAS GPAEGRVARG SERFRDLVPN YNPDIIFKDE ENSGADRLMT ERCKERVNAL AIAVMNMWPG VRLRVTEGWD EDGHHAQDSL HYEGRALDIT TSDRDRNKYG LLARLAVEAG FDWVYYESRN HVHVSVKADN

image



Background

Alternative Names	Desert hedgehog protein (DHH) (HHG-3) [Cleaved into: Desert hedgehog protein N-product; Desert hedgehog protein C-product]
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human DHH Protein
Catalog #: PHH4349

