


Product Name: Recombinant Human ITM2A Protein
Catalog #: PHH4352



Summary

Name	ITM2A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	O43736
Host	HEK293
Species	Human
Molecular Mass	63 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KSTIYRGEMC FFDSEDPANS LRGGEPNFLP VTEEADIREDDNIAIIDVPV PSFSDSDPAA IIHDFEKGMT AYLDLLLGNC YLMPLNTSIV MPPKNLVELF GKLASGRYLP QTYVVREDLV AVEEIRDVSN LGIFIYQLCN NRKSFRLRRR DLLLGFNKRA IDKCWKIRHF PNEFIVETKI

image



Background

Alternative Names Integral membrane protein 2A (Protein E25)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ITM2A Protein
Catalog #: PHH4352

