



Product Name: Recombinant Human SLIT3 Protein
Catalog #: PHH4356



Summary

Name	SLIT3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	O75094
Host	HEK293
Species	Human
Molecular Mass	60 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SAASVDCHGL GLRAVPRGIP RNAERLDLDR NNITRITKMD FAGLKNLRVL HLEDNQVSVI ERGAFQDLKQ LERLRNKNK LQVLPELLFQ STPKLTRLDL SENQIQGIPR KAFRGITDVK NLQLDNNHIS CIEDGAFRAL RDLEILTLNN

image



Background

Alternative Names	Slit homolog 3 protein (Slit-3) (Multiple epidermal growth factor-like domains protein 5) (Multiple EGF-like domains protein 5)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use