


Product Name: Recombinant Human FGF19 Protein
Catalog #: PHH4365



Summary

Name	FGF19_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	tag free
Accession #	O95750
Host	HEK293
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LAFSDAGPHV HYGWGDPIRL RHLYTSGPHG LSSCFLRIRA DGVVDCARGQ SAHSLLEIKA VALRTVAIKG VHSVRYLCMG ADGKMQGLLQ YSEEDCAFEE EIRPDGYNVY RSEKHRLPVS LSSAKQRQLY KNRGFLPLSH FLPMLPMVPE EPEDLRGHLE SDMFSSPLET DSMDPFGLVT GLEAVRSPSF EK

image



Background

Alternative Names Fibroblast growth factor 19 (FGF-19)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human FGF19 Protein
Catalog #: PHH4365

