


Summary

Name	FUCO_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P04066
Host	HEK293
Species	Human
Molecular Mass	51 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QPPRRYTPDW PSLDSRPLPA WFDEAKFGVF IHWGVFSVPA WGSEWFWWHW QGEGRPQYQR FMRDNYPPGF SYADFGPQFT ARFFHPEEWA DLFQAAGAKY VVLTTKHHEG FTNWPSVSW NWNSKDVGPH RDLVGELGTA LRKRNIRYGL YHSLEWFHP LYLLDKKNGF KTQHFVSAKT MPELYDLVNS YKPDLIWSDG EWECPTYWN STNFLSWLYN DSPVKDEVVV NDRWGQNCSC HHGGYYNCED KFKPQSLPDH KWEMCTSIDK FSWGYYRRDMA LSDVTEESI ISELVQTVSL GGNYLLNIGP TKDGLVPIF QERLLAVGKW LSINGEAIYA SKPWRVQWEK NTTSVWYTSK GSAVYAIFLH WPENGLVNL SPITSTTKI TMLGIQGLK WSTDPDKGLF ISLPQLPPSA VPAEFAWTIK LTGVK

image



Background

Product Name: Recombinant Human FUCA1 Protein
Catalog #: PHH4377



Alternative Names

Tissue alpha-L-fucosidase (EC 3.2.1.51) (Alpha-L-fucosidase I) (Alpha-L-fucoside fucohydrolase 1) (Alpha-L-fucosidase 1)

Note

For Research Use Only , Not for Diagnostic Use