

Product Name: Recombinant Human MMP10 Protein
Catalog #: PHH4395



Summary

Name	MMP10_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P09238
Host	HEK293
Species	Human
Molecular Mass	53 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YPLSGAAKEE DSNKDLAQQY LEKYYNLEKD VKQFRRKDSN LIVKKIQGMQ KFLGLEVTGK LDTDTLEVMR KPRCGVPDVG HFSSFPGMPK WRKTHLTYRI VNYTPDLPRD AVDSAIEKAL KVVEEVTPLT FSRLYEGEAD IMISFAVKEH GDFYSFDGPG HSLAHAYPPG PGLYGDIHFD DDEKWTEDES GTNLFLVAAH ELGHSLGLFH SANTEALMYP LYNSFTELAQ FRLSQDDVNG IQSLYGPPPA STEEPLVPTK SVPSGSEMPA KCDPALSFDA ISTLRGEYLF FKDRYFWRRS HWNPEPEFHL ISAFWPSLPS YLDAAYEVNS RDTVFIKGN EFWAIRGNEV QAGYPRGIHT LGFPPTIRKI DAAVSDKEKK KTYFFAADKY WRFDENSQSM EQGFPRLIAD DFPGVPEPKVD AVLQAFGFFY FFSGSSQFEF DPNARMVTHI LKSNSWLHC

image



Background

Product Name: Recombinant Human MMP10 Protein
Catalog #: PHH4395



Alternative Names

Stromelysin-2 (SL-2) (EC 3.4.24.22) (Matrix metalloproteinase-10) (MMP-10)
(Transin-2)

Note

For Research Use Only , Not for Diagnostic Use