



Product Name: Recombinant Human OSM Protein
Catalog #: PHH4413



Summary

Name	ONCM_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P13725
Host	HEK293
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AAIGSCSKEY RVLLGQLQKQ TDLMQDTSRL LDPYIRIQGL DVPKLRHCR ERPGAFPSEE TLRGLGRRGF LQTLNATLGC VLHRLADLEQ RLPKAQDLER SGLNIEDLEK LQMARNILG LRNNIYCMAQ LLDNSDTAEP TKAGRGASQP PTTPASDAF QRKLEGCRFL HGYHRFMHSV GRVFSKWGES PNRSRR

image



Background

Alternative Names Oncostatin-M (OSM)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human OSM Protein
Catalog #: PHH4413

