


**Product Name: Recombinant Human BGN Protein**  
**Catalog #: PHH4430**



## Summary

<b>Name</b>	PGS1_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P21810
<b>Host</b>	HEK293
<b>Species</b>	Human
<b>Molecular Mass</b>	38 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	DEEASGADTS GVLDPDSVTP TYSAMCPFGC HCHLRVVQCS DLGLKSVPE ISPDTLLDL QNNDISELRK DDFKGLQHLY ALVLVNNKIS KIHEKAFSPL RKLQKLYISK NHLVEIPPNL PSSLVELRIH DNRIRKVPKG VFSGLRNMNC IEMGGNPLEN SGFEPGAFDG LKLNLYRISE AKLTGIPKDL PETLNELHLD HNKIQAI ELE DLLRYSKLYR LGLGHNQIRM IENGSLFLP TLRELHLDNN KLARVPSGLP DLKLLQVVYL HSNNITKGVV NDFCPMGFGV KRAYNGISL FNNPVPYWEV QPATFRCVTD RLAIQFGNYK K

## image



## Background

**Alternative Names** Biglycan (Bone/cartilage proteoglycan I) (PG-S1)

**Product Name: Recombinant Human BGN Protein**  
**Catalog #: PHH4430**

---



**Note**

For Research Use Only , Not for Diagnostic Use