



Product Name: Recombinant Human BGN Protein (Biotin)
Catalog #: PHH4431



Summary

Name	PGS1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P21810
Host	HEK293
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DEEASGADTS GVLDPDSVTP TYSAMCPFGC HCHLRVVQCS DLGLKSVPKE ISPDTTLLDL QNNDISELRK DDFKGLQHLY ALVLVNNKIS KIHEKAFSPL RKLQKLYISK NHLVEIPPNL PSSLVELRIH DNRIRKVPKG VFSGLRNMNC IEMGGNPLEN SGFEPGAFDG LKLNLYRISE AKLTGIPKDL PETLNELHLD HNKIQAIELE DLLRYSKLYR LGLGHNQIRM IENGSLSLFP TLRELHLDNN KLARVPSGLP DLKLLQVVYL HSNNITKVG V ND FCPMGFGV KRAYNGISL FNNPVPYWEV QPATFRCVTD RLAIQFGNYK K

image



Background

Alternative Names Biglycan (Bone/cartilage proteoglycan I) (PG-S1)



Product Name: Recombinant Human BGN Protein (Biotin)
Catalog #: PHH4431

Note

For Research Use Only , Not for Diagnostic Use