



Summary

| | |
|------------------------------|--|
| Name | PTPRD_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | P23468 |
| Host | HEK293 |
| Species | Human |
| Molecular Mass | 139 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | VDQTGVSGGV ASFICQATGD PRPKIVWNKK GKKVSNQRFE VIEFDDGSGS VLRIQPLRTP RDEAIYECVA SNNVGEISVS TRLTVLREDQ IPRGFPTIDM GPQLKVVERT RTATMLCAAS GNPDPPEITWF KDFLPVDTSN NNGRIKQLRS ESIGGTPIRG ALQIEQSEES DQGKYECVAT NSAGTRYSAP ANLYVRELRE VRRVPPRFSI PPTNHEIMPG GSVNITCVAV GSPMPYVKWM LGAEDLTPED DMPIGRNVLE LNDVRQSANY TCVAMSTLGV IEAIAQITVK ALPKPPGTPV VTESTATSIT LTWDSGNPEP VSYIIQHQP KNSEELYKEI DGVATTRYSV AGLSPYSYDYE FRVVAVNNIG RGPPSEPVLT QTSEQAPSSA PRDVQARMLS STTILVQWKE PEEPNGQIQG YRVYYTMDPT QHVNNWMKHN VADSQITTIG NLVPQKTYSV KVLAFSIGD GPLSSDIQVI TQTGVPGQPL NFKAEPES SILLSWTPPR SDTIANYELV YKDGEHGEEQ RITIEPGTSY RLQGLKPNLS YFRLAARSP QGLGASTAEI SARTMQSKPS APPQDISCTS PSSTSILVSW QPPPVEKQNG IITEYSIKYT AVDGEDDKPH EILGIPSDTT KYLLEQLEKW TEYRITVTAH TDVGPPESEL SVLIRTNEDV PSGPPRKVEV EAVNSTSVKV SWRSPVFNKQ HGQIRGYQVH YVRMENGEPK GQPMLKDVML ADAQWEFDDT TEHDMIISGL QPETSYSLV TAYTTKGDGA RSKPKLVSTT GAVPGKPRLV INHTQMNTAL IQWHPPVDTF GPLQGYRLKF GRKDMEPLTT LEFSEKEDHF TATDIHKGAS YVFRLSARNK VGFGEEMVKE ISIPEEVPTG FPQNLHSEGT TSTSVQLSWQ PPVLAERNGI ITKYTLLYRD INIPLPMEQ LIVPADTTMT LTGLKPDITY DVKVRHTSK GPGYSPSVQ FRTLTPVDQVF AKNFHVKAVM KTSVLLSWEI PENYNSAMPF KILYDDGKMV EEVDGRATQK LIVNLKPEKS YSFVLTNRGN SAGGLQHRVT AKTAPDVLRT KPAFIGKTNL DGMITVQLPE VPANENIKGY YIIIVPLKKS RGKFIKPWES PDEMEDELLE KEISRKRSSI RYGREVELKP YIAAHFDVLP TEFTLGDDKH YGGFTNKQLQ SGQEYVFFVL AVMEHAESKM YATSPYSDPV VSMDLDPQPI |

Product Name: Recombinant Human PTPRD Protein
Catalog #: PHH4440



image



Background

Alternative Names

Receptor-type tyrosine-protein phosphatase delta (Protein-tyrosine phosphatase delta) (R-PTP-delta) (EC 3.1.3.48)

Note

For Research Use Only , Not for Diagnostic Use