


Product Name: Recombinant Human CTSS Protein
Catalog #: PHH4446



Summary

Name	CATS_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P25774
Host	HEK293
Species	Human
Molecular Mass	78 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QLHKDPTLDH HWHLWKKTYG KQYKEKNEEA VRRLIWEKNL KFMLHNLH SMGMHSYDLG MNHLGDMTSE EVMSLMSSLR VPSQWQRNIT YKSNPNRILP DSVDWREKGC VTEVKYQGSC GACWAFSAVG ALEAQLKLT GKLVSLSAQN LVDCSTEKYG NKGCGGGFMT TAFQYIIDNK GIDSDASYPY KAMDQKQCQYD SKYRAATCSK YTELPYGRED VLKEAVANKG PVSQVVDARH PSFFLYRSGV YYEPSCTQNV NHGVLVVGYG DLNGKEYWLV KNSWGHNFGE EGYIRMARNK GNHCGIASFP SYPEI

image



Background

Alternative Names Cathepsin S (EC 3.4.22.27)

Product Name: Recombinant Human CTSS Protein
Catalog #: PHH4446



Note

For Research Use Only , Not for Diagnostic Use