



Product Name: Recombinant Human CBLIF Protein
Catalog #: PHH4456



Summary

Name	IF_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P27352
Host	HEK293
Species	Human
Molecular Mass	44 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	STQTQSSCSV PSAQEPLVNG IQVLMENSVT SSAYPNPSIL IAMNLAGAYN LKAQKLLTYQ LMSSDNNDLT IGQLGLTIMA LTSSCRDPGD KVSILQRQME NWAPSSPNAE ASAFYGPSLA ILALCQKNSE ATLPIAVRFA KTLANSSPF NVDTGAMATL ALTCMYNKIP VGSEEGYRSL FGQVLKDIVE KISMKIKDNG IIGDIYSTGL AMQALSVTPE PSKKEWNCKK TTDMLNEIK QGKFHNPMSI AQILPSLKGK TYLDVPQVTC SPDHEVQPTL PSNPGPGPTS ASNITVIYTI NNQLRGVELL FNETINVSVK SGSVLLVVLE EAQRKNPMFK FETTMTSWGL VVSSINNIAE NVNHKTYWQF LSGVTPLNEG VADYIPFNHE HITANFTQY

image



Background

Alternative Names	Cobalamin binding intrinsic factor (Gastric intrinsic factor) (Intrinsic factor)
--------------------------	--

Product Name: Recombinant Human CBLIF Protein
Catalog #: PHH4456



(IF) (INF)

Note

For Research Use Only , Not for Diagnostic Use