



Summary

Name	GGT5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-9×His
Accession #	P36269
Host	HEK293
Species	Human
Molecular Mass	61 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RHQAPCGPQA FAHAAVAADS KVCSDIGRAI LQQQGSPVDA TIAALVCTSV VNPQSMGLGG GVIFTIYNVT TGKVEVINAR ETVPASHAPS LLDQCAQALP LGTGAQWIGV PGELRGYAEA HRRHGRLPWA QLFQPTIALL RGGHVVPVPL SRFLHNSILR PSLQASTLRQ LFFNGTEPLR PQDPLPWPAL ATTLETVATE GVEVFYTGRL GQMLVEDIAK EGSQTLQDL AKFQPEVVDA LEVPLGDYTL YSPPPPAGGA ILSFILNVLR GFNFSTESMA RPEGRVNVYH HLVELTKFAK GQRWRLGDPR SHPKLQNASR DLLGETLAQL IRQQIDGRGD HQLSHYSLAE AWGHGTGTSH VSVLGEDGSA VAATSTINTP FGAMVYSPT GIILNNELLD LCERCPRGSG TTPSPVSGDR VGGAPGRCWP PVPGERSPSS MVPSILINKA QGSKLVIGGA GGELIISAVA QAIMSKLWLG FDLRAAIAAP ILHVNSKGCV EYEPNFSQEV QRGLQDRGQN QTQRPFLLNV VQAVSQEGAC VYAVSDLRKS

image



Background

Product Name: Recombinant Human GGT5 Protein
Catalog #: PHH4473



Alternative Names

Glutathione hydrolase 5 proenzyme (EC 3.4.19.13) (Gamma-glutamyl transpeptidase-related enzyme) (GGT-rel) (Gamma-glutamyltransferase 5) (GGT 5) (EC 2.3.2.2) (Gamma-glutamyltransferase-like activity 1) (Gamma-glutamyltranspeptidase 5) (Leukotriene-C4 hydrolase) (EC 3.4.19.14) [Cleaved into: Glutathione hydrolase 5 heavy chain; Glutathione hydrolase 5 light chain]

Note

For Research Use Only , Not for Diagnostic Use