



Summary

Name	EPHA5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P54756
Host	HEK293
Species	Human
Molecular Mass	62 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CYSAPRRAPL WTCLLLCAAL RTLLASPSNE VNLLDSRTVM GDLGWIAFPK NGWEEIGEVD ENYAPIHTYQ VCKVMEQNQN NWLLTSWISN EGASRIFIEL KFTLRDCNSL PGGLGTCKET FNMYFESDD QNGRNIKENQ YIKIDTIAAD ESFTELDLGD RVMKLNTEVR DVGPLSKKGF YLAFQDVGAC IALVSVRVYY KKCPSVVRHL AVFPDTITGA DSSQLLEVSG SCVNHSVTDE PPKMHCSAEG EWLVPIGKCM CKAGYEEKNG TCQVCRPGFF KASPHIQSCG KCPPHSYTHE EASTSCVCEK DYFRRES DPP TMACTRPPSA PRNAISNVNE TSVFLEWIPP ADTGGRKDVS YYIACKKNS HAGVCEECGG HVRYLPRQSG LKNTSVMMVD LLAHTNYTFE IEAVNGVSDL SPGARQYVSV NVTTNQAAPS PVTNVKKGKI AKNSISLSWQ EPDRPNGIIL EYEIKYFEKD QETSYTIKS KETTITA EGL KPASVYVFQI RARTAAGYGV FSRRFEFETT PVFAASSDQS

image



Background

Product Name: Recombinant Human EPHA5 Protein
Catalog #: PHH4496



Alternative Names

Ephrin type-A receptor 5 (EC 2.7.10.1) (Brain-specific kinase) (EPH homology kinase 1) (EHK-1) (EPH-like kinase 7) (EK7) (hEK7)

Note

For Research Use Only , Not for Diagnostic Use