



Product Name: Recombinant Human SULT2A1 Protein
Catalog #: PHH4518



Summary

Name	ST2A1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-His
Accession #	Q06520
Host	HEK293
Species	Human
Molecular Mass	35 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSDDFLWFEG IAFPTMGFRS ETLRKVRDEF VIRDEDVIL TYPKSGTNWL AEILCLMHSK GDAKWIQSVP IWERSPWVES EIGYTALSET ESPRLFSSHL PIQLFPKSFF SSKAKVIYLM RNPRDVLVSG YFFWKNMKFI KKPKSWE EYF EWFCQGTVLY GSWFDHIHGW MPMREEKNFL LLSYEELKQD TGRTIEKICQ FLGKTLEPEE LNLILKNSSF QSMKENKMSN YSLLSVDYVV DKAQLLRKGV SGDWKNHFTV AQAEDFDKLF QEKMADLPRE LFPWE

image



Background

Alternative Names	Sulfotransferase 2A1 (ST2A1) (EC 2.8.2.2) (Bile salt sulfotransferase) (EC 2.8.2.14) (Dehydroepiandrosterone sulfotransferase) (DHEA-ST) (DHEA-ST8) (Hydroxysteroid Sulfotransferase) (HST) (ST2) (SULT2A3)
--------------------------	---

Product Name: Recombinant Human SULT2A1 Protein
Catalog #: PHH4518



Note

For Research Use Only , Not for Diagnostic Use