



Summary

Name	EST3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6UWW8
Host	HEK293
Species	Human
Molecular Mass	61 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TGPEVAQPEV DTTLGRVRGR QVGKGTDRL VNVFLGIPFA QPPLGPDFRS APHPAQPWEG VRDASTAPPM CLQDVESMNS SRFVLNGKQQ IFSVSEDCLV LNVYSPAEPV AGSGRPVMVW VHGGALITGA ATSYDGSALA AYGDVVVVTV QYRLGVLGFF STGDEHAPGN QGFLDVVAAL RWVQENIAPF GDLNLCVTVF GGSAGGSIIS GLVLSVAAG LFHRAITQSG VITTPGIIDS HPWPLAQKIA NTLACSSSSP AEMVQCLQQK EGEELVLSKK LKNTIYPLTV DGTVPKSPK ELLKEKPFHS VPFLMGVNNH EFSWLIPRGW GLLDTMEQMS REDMLAISTP VLTSLDVPPE MMPTVIDEYL GSNSDAQAC QAFQEFMGDV FINVPTVSFS RYLRDSGSPV FFYEFQHRPS SFAKIPAWV KADHGAEGAF VFGGPFLMDE SSRLAFPEAT EEEKQLSLTM MAQWTHFART GDPNSKALPP WPQFNQAEQY LEINPVPRAG QKFREAWMQF WSETLPSKIQ QWHQKQKNRK AQEDL

image



Background

Product Name: Recombinant Human CES3 Protein
Catalog #: PHH4559



Alternative Names

Carboxylesterase 3 (EC 3.1.1.1) (Liver carboxylesterase 31 homolog)

Note

For Research Use Only , Not for Diagnostic Use