


Summary

Name	UNC5B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IZJ1
Host	HEK293
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GSEVLPDSFP SAPAEPLPYF LQEPQDAYIV KNKPVELRCR AFPATQIYFK CNGEWVSQND HVTQEGLDEA TGLRVREVQI EVSRQQVEEL FGLEDYWCQC VAWSSAGTTK SRRAYVRIAY LRKNFDQEPL GKEVPLDHEV LLQCRPPEGV PVAEVEWLKN EDVIDPTQDT NFLLTIDHNL IIRQARLSDT ANYTCVAKNI VAKRRSTTAT VIVYVNGGWS SWAEWSPCSN RCGRGWQKRT RTCTNPAPLN GGAFCEGQAF QKTACTIONICP VDGAWTEWSK WSACSTECAH WRSRECMAPP PQNGGRDCSG TLLDSKNCTD GLCMQNKKTL

image



Background

Alternative Names	Netrin receptor UNC5B (Protein unc-5 homolog 2) (Protein unc-5 homolog B) (p53-regulated receptor for death and life protein 1) (p53RDL1)
--------------------------	---

Product Name: Recombinant Human UNC5B Protein
Catalog #: PHH4564



Note

For Research Use Only , Not for Diagnostic Use