



Summary

Name	ERP44_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-9×His
Accession #	Q9BS26
Host	HEK293
Species	Human
Molecular Mass	45 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ITSLDTENID EILNNADVAL VNFYADWCRF SQMLHPIFEE ASDVIKEEFP NENQVVFARV DCDQHSDIAQ RYRISKYPTL KLFRNGMMM K REYRGQRSVK ALADYIRQQK SDPIQEIRDL AEITTLDRSK RNIIGYFEQK DSDNYRVFER VANILHDDCA FLSAFGDVSK PERYSGDNII YKPPGHSAPD MVYLGAMTNF DVTYNWIQDK CVPLVREITF ENGEELTEEG LPFLILFHMK EDTESLEIFQ NEVARQLISE KGTINFLHAD CDKFRHPLLH IQKTPADCPV IAIDFRHMY VFGDFKDVLI PGKLGQFVFD LHSGLHREF HHGPDPTDTA PGEQAQDVAS SPPESSFQKL APSEYRYTLL

image



Background

Alternative Names Endoplasmic reticulum resident protein 44 (ER protein 44) (ERp44)

Product Name: Recombinant Human ERP44 Protein
Catalog #: PHH4600



(Thioredoxin domain-containing protein 4)

Note

For Research Use Only , Not for Diagnostic Use