


Product Name: Recombinant Human ULBP2 Protein
Catalog #: PHH4604



Summary

Name	ULBP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BZM5
Host	HEK293
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GRADPHSLCY DITVIPKFRP GPRWCAVQGG VDEKTFLLHYD CGNKTVTPVS PLGKKLNVTT AWKAQNPVLR EVVDILTEQL RDIQLENYTP KEPLTLQARM SCEQKAEGHS SGSWQFSFDG QIFLLFDSEK RMWTTVHPGA RKMKEKWEND KVVAMSFHYF SMGDCIGWLE DFLMGMDSTL EPSAGAPLAM S

image



Background

Alternative Names	UL16-binding protein 2 (ALCAN-alpha) (NKG2D ligand 2) (N2DL-2) (NKG2DL2) (Retinoic acid early transcript 1H)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ULBP2 Protein
Catalog #: PHH4604

