



Product Name: Recombinant Human IL21R Protein
Catalog #: PHH4612



Summary

Name	IL21R_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9HBE5
Host	HEK293
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PDLVCYTDYL QTVICILEMW NLHPSTLTLT WQDQYEELKD EATSCSLHRS AHNATHATYT CHMDVFHFMA DDIFSVNITD QSGNYSQECG SFLAESIKP APPFNVTVTF SGQYNISWRS DYEDPAFYML KGKLQYELQY RNRGDPWAVS PRRKLISVDS RSVSLLPLEF RKDSSYELQV RAGPMPGSSY QGTWSEWSDP VIFQTQSEEL

image



Background

Alternative Names	Interleukin-21 receptor (IL-21 receptor) (IL-21R) (Novel interleukin receptor) (CD antigen CD360)
--------------------------	---

Note

Product Name: Recombinant Human IL21R Protein
Catalog #: PHH4612



For Research Use Only , Not for Diagnostic Use