



## Summary

<b>Name</b>	IRPL1_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q9NZN1
<b>Host</b>	HEK293
<b>Species</b>	Human
<b>Molecular Mass</b>	39 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	VVTKRGSADG CTDWSIDIKK YQVLVGEPVR IKCALFYGYI RTNYSLAQSA GLSLMWYKSS GPGDFEPIA FDGSRMSKEE DSIWFRPTLL QDSGLYACVI RNSTYCMKVS ISLTVGENDT GLCYNSKMKY FEKAELSKSK EISCRDIEDF LLPTREPEIL WYKECRTKTW RPSIVFKRDT LLIREVREDD IGNYTCELKY GGFVVRTTE LVTAPLTDK PPKLLYPMES KLIQETQLG DSANLTCRAF FGYSGDV SPL IYWMKGEKFI EDLDENRVWE SDIRILKEHL GEQEVSLI VDSVEEGDLG NYSCYVENGN GRRHASVLLH

## image



## Background

<b>Alternative Names</b>	Interleukin-1 receptor accessory protein-like 1 (IL-1-RAPL-1) (IL-1RAPL-1) (IL1RAPL-1) (EC 3.2.2.6) (Oligophrenin-4) (Three immunoglobulin domain-
--------------------------	--

**Product Name: Recombinant Human IL1RAPL1 Protein**  
**Catalog #: PHH4627**



---

containing IL-1 receptor-related 2) (TIGIRR-2) (X-linked interleukin-1 receptor accessory protein-like 1)

**Note**

For Research Use Only , Not for Diagnostic Use