



**Product Name: Recombinant Human C1QTNF9B Protein**  
**Catalog #: PHH4687**



## Summary

<b>Name</b>	C1T9B_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	N-MBP & C-His
<b>Accession #</b>	B2RNN3
<b>Host</b>	HEK293
<b>Species</b>	Human
<b>Molecular Mass</b>	58 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	DTCRQGHPI PGNPGHNL P GRDGRDGAKG DKG DAGEPGC PGSPGKDGTS GEKGERGADG KVEAKGIKGD QGSRGSPGKH GPKGLAGPMG EKGLRGETGP QGQKGNKGDV GPTGPEGPRG NIGPLGPTGL PGPMGPIGKP GPKG EAGPTG PQGEPGVRGI

## image



## Background

<b>Alternative Names</b>	Complement C1q and tumor necrosis factor-related protein 9B (C1q/TNF-related protein 9B) (CTRP9B) (Complement C1q and tumor necrosis factor-related protein 9-like)
--------------------------	---

## Note

**Product Name: Recombinant Human C1QTNF9B Protein**  
**Catalog #: PHH4687**



---

For Research Use Only , Not for Diagnostic Use