


Product Name: Recombinant Human PLXNB1 Protein
Catalog #: PHH4721



Summary

Name	PLXB1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	O43157
Host	HEK293
Species	Human
Molecular Mass	43 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SPSPPTAPKA LATPAPDTLP VEPGAPSTAT ASDISPGASP SLLSPWGPWA GSGSISSPGS TGSPLHEEPS PPSPQNGPGT AVPAPTDFRP SATPEDLLAS PLSPSEVAHV PPADPGPEAL HPTVPLDLPP ATPATTFFPG AMGSVKPALD

image



Background

Alternative Names Plexin-B1 (Semaphorin receptor SEP)

Note

For Research Use Only , Not for Diagnostic Use