



Product Name: Recombinant Human FZD6 Protein
Catalog #: PHH4735



Summary

Name	FZD6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O60353
Host	HEK293
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IFLPLLRGHS LFTCEPITVP RCMKMAYNMT FFPNLMGHYD QSIAAVEMEH FLPLANLECS PNIETFLCKA FVPTCIEQIH VVPPCRKLCE KVYSDCKKLI DTFGIRWPEE LECDRQLQYCD ETVPVTFDPH TEFLGPQKKT

image



Background

Alternative Names Frizzled-6 (Fz-6) (hFz6)

Note

For Research Use Only , Not for Diagnostic Use