


**Product Name: Recombinant Human LILRB3 Protein**  
**Catalog #: PHH4740**



## Summary

<b>Name</b>	LIRB3_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	O75022
<b>Host</b>	HEK293
<b>Species</b>	Human
<b>Molecular Mass</b>	47 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	LWAEPGSVIS WGSPVTIWCQ GSQEAQEYRL HKEGSPEPLD RNNPLEPKNK ARFSIPSMTE HHAGRYRCHY YSSAGWSEPS DPLEMVMGTGA YSKPTLSALP SPVVASGGNM TLRCGSQKGY HHFVLMKEGE HQLPRTLDSQ QLHSRGFQAL FPVGPVTPSH RWRFTCYYYY TNTPWVWSHP SDPLEILPSG VSRKPSLLTL QGPVLAPGQS LTLQCGSDVG YNRFVLYKEG ERDFLQRPQG QPQAGLSQAN FTLGPVSPSN GGQYRCYGAH NLSSEWSAPS DPLNILMAGQ IYDTVLSAQ PGPTVASGEN VTLLCQSWWQ FDTFLLTKEG AAHPPLRLRS MYGAHKYQAE FPMSPV TSAH AGTYRCYGSY SSNPHELLSH SEPLELVVSG HSGGSSLPPT GPPSTPGLGR

## image



## Background

**Product Name: Recombinant Human LILRB3 Protein**  
**Catalog #: PHH4740**



---

**Alternative Names**

Leukocyte immunoglobulin-like receptor subfamily B member 3 (LIR-3)  
(Leukocyte immunoglobulin-like receptor 3) (CD85 antigen-like family  
member A) (Immunoglobulin-like transcript 5) (ILT-5) (Monocyte inhibitory  
receptor HL9) (CD antigen CD85a)

**Note**

For Research Use Only , Not for Diagnostic Use