


Product Name: Recombinant Human CPD Protein
Catalog #: PHH4757



Summary

Name	CBPD_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	O75976
Host	HEK293
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VIVPSLNPDG RERAQEKDCT SKIGQTNARG KDLDTDFTNN ASQPETKAI ENLIQKQDFS LVALDGGSM LVTYPYDKPV

image



Background

Alternative Names Carboxypeptidase D (EC 3.4.17.22) (Metalloprotease D) (gp180)

Note

For Research Use Only , Not for Diagnostic Use